



REMOVABLE APPLIANCE INFORMATION (ALL RESIN PARTIAL DENTURE)

To keep your new partial looking and feeling like new, please follow these simple directions:

1. With your new partial we have included 5 FREE POST OPS WITHIN THE FIRST YEAR
2. **Each time you place the appliance in your mouth you MUST run under hot tap water for approximately 10-15 seconds or until it becomes flexible. DO NOT bite appliance in place as it is very fragile.
3. **When not in your mouth, DO NOT store or soak your appliance in water. It should stay dry and remain in a safe place.
4. **AVOID using any denture/partial cleaner, bleach, peroxide, or similar cleaning products for your appliance, for these may affect the natural hue of your partial.
5. Clean your appliance regularly. Simply use dish washing liquid and water on your partial. Rinse well, and let air dry.
6. If possible, after eating, rinse your appliance to remove any food particles.
7. Brushing the appliance is not necessary. However, if food particles remain on your appliance after rinsing, they may be removed with a VERY soft brush.
8. Remember to brush your natural teeth and/or gums regularly as directed.
9. Toothpaste and harsh toothbrushes may scratch or dull the surface.
10. Please keep your appliance out of reach of pets. Our little friends find dental appliances to be attractive chew toys.
11. Your flexible appliance is designed to give long term performance under NORMAL usage. Damage or breakage can occur, so handle with care.

**Please call 317-867-5511 if you experience any discomfort from the appliance.
We urge you to maintain your natural teeth as long as possible.**

Please visit us regularly to monitor and professionally clean your appliance.